

Code	Description	Box Qty
710PUR	2x 10Ch needles and PUR drains with Y connector and bag with tap	25
714PUR	2x 14Ch needles and PUR drains with Y connector and bag with tap	25
691PUR	10Ch PUR drain and bag with tap	25
676PUR	14Ch PUR drain and bag with tap	25
694PUR	18Ch PUR drain and bag with tap	25
591PUR	10Ch Double length PUR drain and bag with tap	25
692PUR	10Ch PUR drain and tapless bag	25
677PUR	14Ch PUR drain and tapless bag	25
695PUR	18Ch PUR drain and tapless bag	25
691	10Ch Silicone drain and bag with tap	25
676	14Ch Silicone drain and bag with tap	25
694	18Ch Silicone drain and bag with tap	25
591	10Ch Double length Silicone drain and bag with tap	25
692	10Ch Silicone drain and tapless bag	25
677	14Ch Silicone drain and tapless bag	25
695	18Ch Silicone drain and tapless bag	25
445	700ml Replacement collection bag with tap	25
446	700ml Replacement collection tapless bag	25

Distributed by:

Summit Medical Limited
 Industrial Park
 Bourton on the Water
 Gloucestershire
 GL54 2HQ
 United Kingdom

Tel: +44 (0) 1451 821311
 Fax: +44 (0) 1451 821092
 Customer Services Tel: +44 (0) 1451 810900
 Customer Services Fax: +44 (0) 1451 822302
 e-mail: info@summit-medical.co.uk
 website: www.summit-medical.co.uk

Summit Medical Operations
 BSI Certificate No. FM 32336
 ISO9001/EN46001



LowVac[®] Plus Wound Drainage System

Effective, Progressive, Innovative



SUMMIT MEDICAL



LowVac[®] Plus Wound Drainage System

The purpose of wound drainage is to promote faster wound healing. Effective wound drainage helps to reduce the incidence of haematomas forming and reduce the possibilities of infection by the removal of air and fluids from the wound site. Wound drainage also helps to release internal pressure and reduce the formation of scar tissue.

LowVac[®] Plus Wound Drainage system

LowVac[®] Plus is a constant low vacuum wound drainage system which can be used in a number of surgical procedures such as:

- Orthopaedic Surgery e.g. knee, hip, spinal, hand
- General Surgery e.g. Abdominal
- Gynaecology
- Plastic Surgery
- Breast Surgery



- **Constant low vacuum drainage**
This system operates below 100mmHg and provides a vacuum during its operation which reduces the incidence of haematomas forming and minimizes the risk of patient infection. It also eliminates the need for constant monitoring of the device.
- **Easy To Use**
Simple one-handed priming bellows and pinch clamps.
- **Different types of drains available**
Silicone Drains - soft to minimise patient trauma.
Polyurethane (PUR) Drains – more haemocompatible than PVC drains, ease of removal (1)
- **Full range of drain sizes available – sizes 10, 14 & 18Ch**
Aids suitability for all types of procedures.
- **Available with Y Connector**
Size 10Ch and 14Ch available. Two pre-assembled trocars and drains included with the device, Y connector already attached.
- **Two non return valves**
Prevents any back flow of exudate into the wound.
Minimises the potential for retrograde infection of the wound.
- **Safety caps on the bag and bellows**
Minimises accidental cross infection and spillage.
- **Space for recording patient details on collection bag**
Allows for clear labelling, saving nursing time.
- **Measuring scale on bag**
Allows assessment of the volume of fluid collected.
- **Drain connectors pre-assembled**
No additional assembly required, saving user time.
- **Effective, easy to use pinch clamps**
System can be primed in line with user preference.
- **Latex free product**
Provides protection against potential latex allergy for nursing staff and patients.
- **PVC free packaging**
Helps minimise environmental pollution.
- **Tap or Tapless 700ml Bags**
Replacement bags available.