

Silicone Reservoir and Closed Systems

FEATURES

Post-operative drainage system with low vacuum aspiration level, ideal for use with Jackson-Pratt type silicone drains. Clear silicone reservoirs available in various capacities and featuring anti-reflux valve. Drentech® Pump standard and Drentech® Closed with collection bag 200 & 400 ml reservoirs are with double connector for allowing multiple drains attachment.

Drentech Pump

CODE	CAPACITY
10400	100 ml
10410	200 ml
10420	400 ml
PACKAGING: 10 PIECES PER BOX	

Drentech Closed

CODE	CAPACITY	BAG CAPACITY
22610	100 ml	200 ml
22620	200 ml	600 ml
22630	200 ml	200 ml
22600	400 ml	600 ml
PACKAGING: 10 PIECES PER BOX		

Replacement bag is available separately from the set, sterile package

Collection bag

CODE	CAPACITY
12004	600 ml
12014	600 ml no tap
PACKAGING: 25 PIECES PER BOX	
12015	200 ml
PACKAGING: 50 PIECES PER BOX	



REDAX®

dall'esperienza, le nuove tecnologie del drenaggio

Via Volta, 51 - 41037 Mirandola (MO) ITALY
tel. +39 0535 611673 fax +39 0535 23224
e-mail: info@redax.it - www.redax.it



Silicone Drainage System



REDAX®

dall'esperienza, le nuove tecnologie del drenaggio


Redax manufactures a wide range of high quality drains designed for post-operative applications.

Silicone drains are developed in different profile and fenestration to satisfy any clinical needs.

All silicone devices are tested for haemocompatibility and biocompatibility.

Drain tubes are designed with atraumatic profile and feature radiopaque x-ray line.

Redax realizes also a complete range of silicone reservoirs for post-operative low vacuum drainage. Silicone bulb are available also with collection bag.




Pre cut Silicone Wound Drains

FEATURES

High quality silicone drain designed mainly for gravity • Pre-cut 50 cm length • Atraumatic elliptical eyes • Radiopaque line along entire length and depth reference marks every 2 cm. Available in a wide range of size also with 2 liter collection bag.

CODE STANDARD	CODE DRAIN WITH BAG	EXT. DIAMETER	
		CH	mm
24309	----	9	3
24312	21712	12	4
24315	21715	15	5
24318	21718	18	6
24321	21721	21	7
24324	21724	24	8
24328	21728	28	9,3
24330	21730	30	10
24333	21733	33	11
24336	----	36	12

20 PIECES PACKAGING: PER BOX
25 PIECES




Penrose Silicone Drains

FEATURES

Capillary silicone tube available in various sizes, length 30 and 50 cm • Radiopaque line along entire length. Double package, sterile.

CODE 30 cm	CODE 50 cm	EXT. DIAMETER mm
22206	22006	6
22208	22008	8
22210	22010	10
22212	22012	12
22214	22014	14
22216	22016	16
22240	22040	40 Open

PACKAGING: 20 PIECES PER BOX




T-Kehr biliary drains

FEATURES

T-shape silicone tubes for biliary drainage. Total length is 50 cm and horizontal part 15 cm length.

CODE STANDARD	EXT. DIAMETER	
	CH	mm
24409	9	3
24412	12	4
24415	15	5
24418	18	6
24421	21	7
24424	24	8

PACKAGING: 10 PIECES PER BOX



Flat and round Silicone Drains Jackson-Pratt® type

FEATURES

Silicone drains like Drains Jackson-Pratt type are realized in two different profile: flat and round with a special fenestration of 20 cm length. Tube is totally transparent, ideal for post-operative drainage allowing continuous aspiration. Made of soft material ensuring atraumatic patient drainage.

FLAT DRAINS		
TOTAL FENESTRATION	FENESTRATION 3/4	SIZE mm
24110	24100	3X7
24111	24101	3x7 w/trocar
24130	24120	4x10
24131	24121	4x10 w/trocar
24140		4X13
24141		4x13 w/trocar


PACKAGING: 20 PIECES PER BOX

ROUND DRAINS		
STANDARD	W/TROCAR	SIZE
24200	24201	2,3 mm - CH 7
24210	24211	3,2 mm - CH 10
24220	24221	4,9 mm - CH 15
24230	24231	6,3 mm - CH 19

PACKAGING: 20 PIECES PER BOX



Inner ribbing in flat tube reduces risk of clogging and collapsing.



Drentech® Bag - Robinson System type

FEATURES

Closed drainage system designed for gravity, with silicone tube. Available with 100 cm silicone drain attached to a collection bag. Silicone tube is radiopaque with atraumatic elliptical eyes. Collection bag has a capacity of 600 ml, equipped with internal one-way valve to prevent reflux.

CODE	EXTERNAL DIAMETER	
	CH	mm
21409	9	3
21412	12	4
21415	15	5
21418	18	6
21420	21	7
21424	24	8
21426	26	8,7
21428	28	9,3
21430	30	10
21433	33	11
21436	36	12

PACKAGING: 10 PIECES PER BOX

