

AquaNASE®

AquaNASE® cannula for HFHOT can be supplied with AquaVENT® breathing systems in a complete Flowkit.

Patients can eat, drink, communicate and mobilise pulmonary secretions without disrupting therapy support.

Heated humidification helps prevent drying of upper airway mucosa increasing patient comfort during oxygen therapy.

Flows closer to a patient's inspiratory demand reduces ambient air entrainment and oxygen dilution.

Use AquaNASE® during patients non-invasive ventilation rest periods for:

- oral hygiene
- hydration
- feeding
- physiotherapy
- managed weaning

| Clinical applications ⁽¹⁾

- Acute hypoxemic respiratory failure
- Post-extubation period
- Pre-intubation
- Emergency department
- Bronchoscopy and other invasive procedures
- Palliative care
- Acute heart failure
- Chronic airway disease

