


## NIV Face Masks for CPAP & BPAP



Code	Description	Dead Space
HSFM10211Y	Adult Large non-vented AAV yellow elbow	205ml
HSFM10214Y	Adult Medium non-vented AAV yellow elbow	195ml
HSFM10217Y	Adult Small non-vented AAV yellow elbow	165ml
HSFM10211	Adult Large non-vented, blue elbow	205ml
HSFM10214	Adult Medium non-vented, blue elbow	195ml
HSFM10217	Adult Small non-vented, blue elbow	165ml

## Spirale® Drug Delivery System

Nebulisation with COVID-19 patients increases the risk of aerosolisation and transmission of infection to healthcare workers. Use of metered-dose inhalers or spacers is preferred where possible.




## The NHS Guidelines for NIV consumables advise:<sup>1</sup>

- Any patient on acute NIV should be managed with a non-vented mask & an exhalation port in the circuit.
- A viral/bacterial filter should be placed in the circuit between the mask & the exhalation port.


Armstrong Medical NIV Breathing Systems are available in heated and non-heated specifications.

## CPAP Systems



Code	Description	Case Qty	Pallet Qty
<b>a.</b>			
AMCP1409-003	<ul style="list-style-type: none"> <li>• Adult Flowkit AquaVENT® heated humidified oxygen circuit</li> <li>• 1.5m heated limb</li> <li>• Humidification limb 0.6m</li> <li>• Humidification chamber</li> </ul>	20	240
<b>b.</b>			
AMCPUK01240	<ul style="list-style-type: none"> <li>• Adult CPAP System (part of)</li> <li>• Y-piece attachment</li> <li>• 5cm patient PEEP valve</li> <li>• Safety valve assembly</li> <li>• Bacterial viral breathing filter</li> <li>• Manometer line assembly (4mm ID)</li> </ul>	20	240
AMCPUK01258	<ul style="list-style-type: none"> <li>• Adult CPAP System (part of)</li> <li>• Y-piece attachment</li> <li>• 5cm patient PEEP valve</li> <li>• Safety valve assembly</li> <li>• Bacterial viral breathing filter</li> <li>• Manometer line assembly (3mm ID for use with FD140 &amp; FD140i)</li> </ul>	40	480

## BPAP Systems



Code	Description	Case Qty	Pallet Qty
<b>a.</b>			
AMVC1792-032	<ul style="list-style-type: none"> <li>• Heated</li> <li>• Adult ventilator circuit</li> <li>• Single limb</li> <li>• Length 1.5m</li> <li>• Exhalation valve with tubing</li> <li>• 0.2µm hydrophobic gas filter</li> <li>• Humidification limb length 0.9m</li> <li>• Auto-fill humidification chamber</li> <li>• Bacterial viral breathing filter</li> </ul>	20	240
<b>b.</b>			
AMVC1792-010	<ul style="list-style-type: none"> <li>• Non-heated</li> <li>• Adult ventilator circuit</li> <li>• Single limb</li> <li>• Length 1.8m</li> <li>• Exhalation valve</li> <li>• Pressure line</li> <li>• Bacterial viral breathing filter</li> </ul>	40	480
AMVC1792-013	<ul style="list-style-type: none"> <li>• Non-heated</li> <li>• Adult ventilator circuit</li> <li>• Single limb</li> <li>• Length 1.8m</li> <li>• Exhalation valve</li> <li>• Bacterial viral breathing filter</li> </ul>	40	480

References:

<sup>1</sup> <https://www.england.nhs.uk/coronavirus/wp-content/uploads/sites/52/2020/03/specialty-guide-NIV-respiratory-support-and-coronavirus-v3.pdf>